

## ■ Announcements

- ❖ Discuss mid-term feedback

- ❖ FNT heads up: methods, 3Q (got rid of one)

- ❖ Day 7: quiz; staggered arrivals (~1 – 1.5hrs)

- ❖ Office hours

## ■ Pre-lab Lecture

- ❖ SDS-PAGE

- ❖ Affinity purification recap

- ❖ Today in Lab (Mod 2 Day 6)

# SDS-PAGE preparation

↖ Acrylamide - toxic


- You will make whole cell extracts with equal cell #s

- Based on OD<sub>600</sub> reading, normalize (1) OD = 1.0 (2) OD = 0.5

V<sub>max</sub> = 15mL

(1) 7.5mL + 7.5mL H<sub>2</sub>O

(2) 15mL



- Gel separates proteins based on size, shape, charge

- Sample preparation

- SDS: coat proteins with negative charge

in hood { – β-Me: breaks S-S bonds


– Boiling: denature higher-order structures

– Sample Buffer has SDS, β-Me, plus glycerol, BPB dye

↑  
make uniform

# SDS-PAGE visualization, analysis

- Visualization: Coomassie stain (binds certain AA)
- Two ladders: visualization, quantification





# Affinity purification

- Basis:

His-tag in vector

6x, binds to metals




Courtesy of Life Technologies, Carlsbad, CA. Used with permission.

high [imidazole]

elute His-tagged protein

# Today in Lab

- Lyse cell pellets in BPER
  - BSA “carrier,” protease inhibitors
  - Add 4 mL lysis enzymes
- Run a 25  $\mu$ L aliquot through SDS-PAGE
  - Two ladders also  $\rightarrow$  boil these too
  - Stick with equal volumes if you have  $< 25$
- Purify IPC protein from the rest (long!)
  - Immediately take 10  $\mu$ L aliquot and measure concentration
  - The rest is stabilized w/BSA, to be titrated against calcium next time



MIT OpenCourseWare  
<http://ocw.mit.edu>

20.109 Laboratory Fundamentals in Biological Engineering  
Spring 2010

For information about citing these materials or our Terms of Use, visit: <http://ocw.mit.edu/terms>.