

PBC Day 6 Interpretation Questions (with answers)

1. Interpret your Western blot data. Your interpretation should include answers to the following questions:

- a. Which lanes have bands and which do not?
answers will vary from student to student
- b. Is this what you expected?
If bands are observed in the lanes mentioned in the answer to part c below, then the answer should be "yes" — otherwise, "no." Many students' blots also had bands at MW smaller than that of full length Bgal (~112 kDa); these are probably degraded fragments of Bgal, and they should mention these if they saw them.
- c. If you don't see any bands, where would you have expected to see them?
Ideally, you would expect to see bands in the following lanes: CL; CL-S; AS-P; DEAE/AF Load; DEAE 0.4M fraction; AF before PD-10, AF after PD-10. Many students will likely see a band(s) in CL-P, as cell lysis was incomplete. These lanes should correspond to fractions that had Bgal activity in their ONPG assay!
- d. Suggest reason(s) that your "expected" bands might be missing.
Some possibilities include: forgetting to load a sample/didn't have a sample to load; low protein amounts in sample (below level of detection of Western blot); poor transfer; or experimental errors (forgot to add NBT/BCIP, etc.)

2. Enzyme Kinetics Analysis for **EACH** strain (CSH and HIS461). There is a guide in the manual appendix (page 62-66) that will help you through these calculations. Remember: everything happens **in the reaction tube**, and DON'T FORGET UNITS!!

Use Excel or some comparable program for data analysis and include all relevant charts and graphs:

- i. Plot of primary data for EACH strain


Will vary, but watch out for whether they remember to multiple A_{420} by 2 (to get the A_{420} in the reaction tube).

PBC Day 6 Interpretation Questions and Answers (continued)

- ii. Michaelis-Menton plot for EACH strain— Estimate K_M and V_{max}


Will vary, but looks like this:

PBC Day 6 Interpretation Questions and Answers (continued)



iii. Lineweaver-Burke plot for EACH strain – Calculate K_M and V_{max}

Will vary, but looks like this:



iv. Calculation of k_{cat} for EACH strain.

Varies but the calculation for $[E_t]$ should look like:

1. $\frac{\text{volume of sample added } (\mu\text{L})}{15 \text{ mL}} \times \frac{1 \text{ mL}}{1000 \mu\text{L}} \times \text{purity}^*(\text{mg/mL}) = [\text{protein}] \text{ in rxn.}$
2. $\text{protein conc. in rxn. (mg/mL)} \times 1 \text{ g}/1000 \text{ mg} = \text{g/mL in rxn.}$
3. $\text{g/mL} \times \text{mol/g (from Bgal MW)} = \text{mol/mL in rxn.}$

*purity comes from Coomassie gel, and should have a value between 0-1)

$k_{cat} = V_{max}/[E_t]$ and has the units U/mol