- Announcements
- Lab Quiz
- Pre-lab Lecture
 - Writing a Figure/Caption
 - In Vitro Transcription
 - Today in Lab: M1D3

Announcements


- No lab next Tuesday (Monday schedule) or Wednesday – see you in a week!
- Next time in lab is packed
 - Very short quiz + pre-lab lecture
- FNT: Lots! Reading and calculations for Day 4, practice figures, writing exercises.

Figures: Style and Scope

Title: concise, informative → gives overall goal/result

- Caption: gives <u>context</u> for result, from big to small
 - Introduce what are we looking at
 - Include just methods needed to understand the result
 - Define all elements (e.g., DNA ladder)
 - Cover primarily <u>facts</u>, not interpretation
 e.g. observed/expected sizes
- ◆ Aesthetics simplicity, clarity → at-a-glance labeling e.g. some ladder band sizes

Figures: Example




result, without over-interpreting

FIGURE 5. CCL21 impacts naive T cell proliferation under conditions of rare Ag-specific T-DC encounters. Cocultures comprising 9% OVA-specific OT-II CD4⁺ T cells, 81% C57BL/6 CD4⁺ T cells, 5% OVA-mDC, and 5% iDC with/without CCL21 were analyzed by flow cytometry at 85 h. A, Sample CFSE histograms are shown for control (left, iDC only) and experimental (right, with OVA-mDC) conditions. B, OT-II cell recovery for all conditions is shown. Average \pm SD for 3 wells per condition. *, bracketed conditions statistically different ($p \le 0.05$). A and B, from one representative of five experiments.


overview of experiment walkthrough of figure

Source: Stachowiak, A., et al. "Homeostatic Lymphoid Chemokines Synergize with Adhesion Ligands to Trigger T and B Lymphocyte Chemokinesis." *J Immunol* 177 (2006): 2340-2348. © American Association of Immunologists, Inc. All rights reserved. This content is excluded from our Creative Commons license. For more information, see http://ocw.mit.edu/fairuse.

SELEX Overview



In vitro transcription



Cartoons from Niles Lecture 2.

PCR vs. IVT

PCR	IVT
DNA plasmid template	DNA fragment (linear)
Primers	not needed!
dNTPs	NTPs
Taq DNA polymerase	T7 RNA Polymerase
Buffer, Mg ions	similar

Today in Lab

- Working with RNA
 - Gloves on, keep clean
- Set up IVT rxns
 - Run for 4 hrs, note your start time up front
 - Stored frozen till next time
 - Return the rest of your DNA, too!
- Meanwhile, journal article discussion
- Also sign up for column conditions

MIT OpenCourseWare http://ocw.mit.edu

20.109 Laboratory Fundamentals in Biological Engineering Spring 2010

For information about citing these materials or our Terms of Use, visit: http://ocw.mit.edu/terms.